

# GARMENTS

## 8 CHAMBER STANDARD LEG SLEEVE



### 8 Chamber Standard Leg Sleeve Size and Part Number

Measurements below are in **Inches** and are the maximum allotted point on limb

Part Number	Size		Circumferences						Lengths	
			Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease
G8-3045-S	Small	Standard	31"	30"	28"	24.5"	21.5"	19.5"	25.5"	27"
G8-3045-M	Medium	Standard	32"	30"	28"	24.5"	21.5"	19.5"	28"	32"
G8-3045-L	Large	Standard	33"	32"	30"	26"	21.5"	19.5"	32"	36"
GN8-3045-S	Small	Narrow	25.5"	25"	22"	20"	18"	19.5"	25.5"	27"
GN8-3045-M	Medium	Narrow	27.5"	25.5"	23"	20"	18"	19.5"	28"	32"
GN8-3045-L	Large	Narrow	28"	26"	24"	21"	18"	19.5"	32"	36"
GW8-3045-S	Small	Wide	35.5"	34"	30.5"	28"	25.5"	21.5"	25.5"	27"
GW8-3045-M	Medium	Wide	36.5"	34"	30"	28"	25.5"	22.5"	28"	32"
GW8-3045-L	Large	Wide	37"	36"	32"	29"	24"	20.5"	32"	36"
GXW8-3045	Medium	Extra Wide	44"	41"	36"	34"	30"	22.5"	29"	32"

Measurements below are in **Centimeters** and are the maximum allotted point on limb

Part Number	Size		Circumferences						Lengths	
			Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease
G8-3045-S	Small	Standard	79	76	71	62	55	49	65	69
G8-3045-M	Medium	Standard	81	76	71	62	55	49	71	81
G8-3045-L	Large	Standard	84	81	76	66	55	49	81	91
GN8-3045-S	Small	Narrow	65	63	56	51	46	49	65	69
GN8-3045-M	Medium	Narrow	70	65	58	51	46	49	71	81
GN8-3045-L	Large	Narrow	71	66	61	53	46	49	81	91
GW8-3045-S	Small	Wide	90	86	77	71	65	55	65	69
GW8-3045-M	Medium	Wide	93	86	76	71	65	57	71	81
GW8-3045-L	Large	Wide	94	91	81	74	61	52	81	91
GXW8-3045	Medium	Extra Wide	112	104	91	86	76	57	74	81

HCPCS: E0667 - Full Leg

