

GARMENTS

4 CHAMBER ADJUSTABLE LEG SLEEVE


4 Chamber Adjustable Leg Sleeve Size and Part Number

Measurements below are in **Inches** and are the maximum allotted point on limb

Part Number	Size		Circumferences						Lengths	
			Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease
GA-3045-H	Half	Adjustable	N/A	N/A	32"	30"	23.5"	22"	19"	NA
GA-3045-S	Small	Adjustable	34"	33"	31"	27.5"	21.5"	19.5"	25.5"	27"
GA-3045-M	Medium	Adjustable	35"	34"	31"	27.5"	21.5"	19.5"	28"	32"
GA-3045-L	Large	Adjustable	37"	36"	33"	29"	21.5"	19.5"	32"	36"
GWA-3045-H	Half	Adj. Wide	N/A	N/A	33.5"	31"	29"	23"	19"	N/A
GWA-3045-S	Small	Adj. Wide	38.5"	37"	33.5"	31"	25.5"	21.5"	25.5"	27"
GWA-3045-M	Medium	Adj. Wide	39.5"	36.5"	33.5"	31"	25.5"	22.5"	28"	32"
GWA-3045-L	Large	Adj. Wide	40"	38"	35"	32"	24"	22.5"	32"	36"
GXWA-3045	Medium	Adj. Extra Wide	47"	44"	39"	37"	30"	22.5"	29"	32"

Measurements below are in **Centimeters** and are the maximum allotted point on limb

Part Number	Size		Circumferences						Lengths	
			Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease
GA-3045-H	Half	Adjustable	N/A	N/A	81	76	60	56	48	NA
GA-3045-S	Small	Adjustable	86	84	78	70	55	50	65	69
GA-3045-M	Medium	Adjustable	89	86	79	70	55	50	71	81
GA-3045-L	Large	Adjustable	94	91	84	74	55	50	81	91
GWA-3045-H	Half	Adj. Wide	N/A	N/A	85	79	74	58	48	NA
GWA-3045-S	Small	Adj. Wide	98	94	85	79	65	55	65	69
GWA-3045-M	Medium	Adj. Wide	100	93	85	79	65	57	71	81
GWA-3045-L	Large	Adj. Wide	102	97	89	81	61	57	81	91
GXWA-3045	Medium	Adj. Extra Wide	119	112	99	94	76	57	74	81



HCPCS:

E0669 - Half Leg

E0667 - Full Leg