

# GARMENTS

## 4 CHAMBER ARM SLEEVE

### 4 Chamber Arm Sleeve Size and Part Number

Measurements below are in **Inches** and are the maximum allotted point on limb

Part Number	Size	Circumferences					Lengths	
		Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips
GS-3035-S	Small	23.5"	22"	19.5"	17.5"	15"	24"	26.5"
GS-3035-M	Medium	27"	26"	24"	21"	18"	28"	31"
GS-3035-L	Large	25.5"	24"	21.5"	19"	16"	31"	35.5"

Measurements below are in **Centimeters** and are the maximum allotted point on limb

Part Number	Size	Circumferences					Lengths	
		Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips
GS-3035-S	Small	60	56	49	44	38	61	67
GS-3035-M	Medium	69	66	61	53	46	71	79
GS-3035-L	Large	65	61	55	48	41	79	90

HCPCS: E0668

