

# GARMENTS

## 4 CHAMBER ARM & SHOULDER SLEEVE

### 4 Chamber Arm & Shoulder Sleeve Size and Part Number

Measurements below are in **Inches** and are the maximum allotted point on limb

Part Number	Size	Circumferences					Lengths	
		Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips
GS-3035-SH-S-L	Small Left	23.5"	22"	19.5"	17.5"	15"	24"	26.5"+8"
GS-3035-SH-S-R	Small Right	23.5"	22"	19.5"	17.5"	15"	24"	26.5"+8"
GS-3035-SH-M-L	Medium Left	27"	26"	24"	21"	18"	28"	31"+8"
GS-3035-SH-M-R	Medium Right	27"	26"	24"	21"	18"	28"	31"+8"
GS-3035-SH-L-L	Large Left	25.5"	24"	21.5"	19"	16"	31"	35.5"+8"
GS-3035-SH-L-R	Large Right	25.5"	24"	21.5"	19"	16"	31"	35.5"+8"

Measurements below are in **Centimeters** and are the maximum allotted point on limb

Part Number	Size	Circumferences					Lengths	
		Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips
GS-3035-SH-S-L	Small Left	60	56	49	44	38	61	67+20
GS-3035-SH-S-R	Small Right	60	56	49	44	38	61	67+20
GS-3035-SH-M-L	Medium Left	69	66	61	53	46	71	79+20
GS-3035-SH-M-R	Medium Right	69	66	61	53	46	71	79+20
GS-3035-SH-L-L	Large Left	65	61	55	48	41	79	90+20
GS-3035-SH-L-R	Large Right	65	61	55	48	41	79	90+20

HCPCS: E0668

