

GARMENTS


16 CHAMBER ELITE BILATERAL BIO VEST



HPCS: E0657 – Chest E0668 – Arm



Part Number: GV-3010-S-2/M-2/L-2



Upper Extremity/16 Chamber Garment Sizing Chart																	
Measurements below are in Inches and are the maximum allotted point on limb										Measurements below are in Centimeters and are the maximum allotted point on limb							
		Circumferences					Lengths			Circumferences					Lengths		
Part Number	Size	Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips	Thorax (Adjustable)	Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips	Thorax (Adjustable)
GV-3010-S-L	Small Left	25"	23.5"	21"	20"	19"	25.5"	29.5"	28" to 55"	63	60	53	51	48	65	75	71 to 140
GV-3010-S-R	Small Right	25"	23.5"	21"	20"	19"	25.5"	29.5"	28" to 55"	63	60	53	51	48	65	75	71 to 140
GV-3010-M-L	Medium Left	25"	23.5"	21"	20"	19"	28"	32"	28" to 55"	63	60	53	51	48	71	81	71 to 140
GV-3010-M-R	Medium Right	25"	23.5"	21"	20"	19"	28"	32"	28" to 55"	63	60	53	51	48	71	81	71 to 140
GV-3010-L-L	Large Left	25"	23.5"	21"	20"	19"	31"	35"	28" to 55"	63	60	53	51	48	79	89	71 to 140
GV-3010-L-R	Large Right	25"	23.5"	21"	20"	19"	31"	35"	28" to 55"	63	60	53	51	48	79	89	71 to 140