

GARMENTS

4 CHAMBER STANDARD LEG SLEEVE



4 Chamber Standard Leg Sleeve Size and Part Number

Measurements below are in **Inches** and are the maximum allotted point on limb

Part Number	Size		Circumferences						Lengths	
			Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease
GN-3045-S	Small	Narrow	25.5"	25"	22"	20"	18"	19.5"	25.5"	27"
GN-3045-M	Medium	Narrow	27.5"	25.5"	23"	20"	18"	19.5"	28"	32"
GN-3045-L	Large	Narrow	28"	26"	24"	21"	18"	19.5"	32"	36"
GS-3045-H	Half	Standard	N/A	N/A	29"	27"	23.5"	22"	19"	NA
GS-3045-S	Small	Standard	31"	30"	28"	24.5"	21.5"	19.5"	25.5"	27"
GS-3045-M	Medium	Standard	32"	30"	28"	24.5"	21.5"	19.5"	28"	32"
GS-3045-L	Large	Standard	33"	32"	30"	26"	21.5"	19.5"	32"	36"
GW-3045-H	Half	Wide	N/A	N/A	30.5"	28"	29"	23"	19"	NA
GW-3045-S	Small	Wide	35.5"	34"	30.5"	28"	25.5"	21.5"	25.5"	27"
GW-3045-M	Medium	Wide	36.5"	34"	30"	28"	25.5"	22.5"	28"	32"
GW-3045-L	Large	Wide	37"	36"	32"	29"	24"	20.5"	32"	36"
GXW-3045	Medium	Extra Wide	44"	41"	36"	34"	30"	22.5"	29"	32"

Measurements below are in **Centimeters** and are the maximum allotted point on limb

Part Number	Size		Circumferences						Lengths	
			Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease
GN-3045-S	Small	Narrow	65	64	56	51	46	50	65	69
GN-3045-M	Medium	Narrow	70	65	58	51	46	50	71	81
GN-3045-L	Large	Narrow	71	66	61	53	46	50	81	91
GS-3045-H	Half	Standard	N/A	N/A	74	69	60	56	48	NA
GS-3045-S	Small	Standard	79	76	71	62	55	50	65	69
GS-3045-M	Medium	Standard	81	76	71	62	55	50	71	81
GS-3045-L	Large	Standard	84	81	76	66	55	50	81	91
GW-3045-H	Half	Wide	N/A	N/A	78	71	74	58	48	NA
GW-3045-S	Small	Wide	90	86	77	71	65	55	65	69
GW-3045-M	Medium	Wide	93	86	76	71	65	57	71	81
GW-3045-L	Large	Wide	94	91	81	74	61	52	81	91
GXW-3045	Medium	Extra Wide	112	104	91	86	76	57	74	81

HCPCS: E0669 – Half leg E0667 – Full leg

